


SMALL VICE GRIP

applicable for packaging materials, paper, plastic foils, wires, yarns, etc.


Item no.:	AD 0005
Max. load:	100 N
Opening:	0-8 mm (depending on jaws)
Coupling:	M6 internal thread Further couplings on request
Body:	Aluminium, anodized
Temperature range:	0 ... +130°C (depending on jaws)
Weight:	60 gram per grip without jaws
Scope of supply:	1 piece of grip without jaws



Jaws for AD 0005: Scope of delivery 1 pair = 2 jaws

Item No.:	Surface:	Clamping surface H x W:	Opening:	Temperature range:	Weight per pair
AD 0006	Blank jaws	10 x 25 mm	0 - 8 mm	-70...+280°C	18 Gram
AD 0007	Rubber jaws	10 x 25 mm	0 - 6 mm	-20...+80°C	19 Gram
AD 0008	Pyramid jaws	10 x 25 mm	0 - 8 mm	-70...+280°C	17 Gram
AD 0009	Blank jaws	10 x 50 mm	0 - 6 mm	-70...+280°C	46 Gram
AD 0010	Rubber jaws	10 x 50 mm	0 - 4 mm	-20...+80°C	48 Gram
AD 0011	Pyramid jaws	10 x 50 mm	0 - 6 mm	-70...+280°C	43 Gram
AD 0012	Blank jaws	10 x 100 mm	0 - 6 mm	-70...+280°C	92 Gram
AD 0013	Rubber jaws	10 x 100 mm	0 - 4 mm	-20...+80°C	95 Gram
AD 0014	Pyramid jaws	10 x 100 mm	0 - 6 mm	-70...+280°C	86 Gram

Jaws with other dimensions and surfaces on request



Blanke Backen
- Suitable for own
further treatment
- Steel-blanks without
coating



Gummi Backen
- 1mm rubber coating (NBR)
- Nickel-plated



Pyramiden Backen
- Pyramids 0.8x45°
- Hardened steel 58 HRC
- Nickel-plated